PL

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EN

Root canals filling and sealing material

**BEFORE USE PLEASE READ THE FOLLOWING
INSTRUCTIONS CAREFULLY
PRODUCT FOR DENTAL USE ONLY****COMPOSITION:** epoxy resin, amino resin, inorganic filler, calcium tungstate, excipients.**INTENDED USE**

Permanent root canal filling of the permanent teeth using gutta-percha points. Permanent root canal filling of the permanent teeth using root canal points. The material can be used for all techniques of the root canal sealing. It is perfect for thermic methods.

INSTRUCTIONS FOR USE

The sealing paste SYNTEX comes in a syringe dispensing two pastes in the ratio 1:1.

Before the first use squeeze a small amount of product directly from the syringe to ensure even pastes dispensing and then place the mixing tip onto the syringe.

Product mixing is possible as well without using any mixing tip. In order to do this, squeeze the pastes in the same ratio directly on the mixing paper and mix using any tool until obtaining the uniform color.

Prepare, rinse and dry the root canal prior to application.

Place a small amount of SYNTEX on the tip of Lentulo needle; slowly move Lentulo needle towards the apex at the lowest possible micro-engine speed. Remove Lentulo from the canal using the minimal engine speed. It is possible to apply SYNTEX paste to a previously disinfected and dried point. Place it in the canal with a pumping movement. Additional points may be placed in the canal or it may be filled with the method of lateral or vertical condensation of gutta-percha.

If pulp chamber is contaminated with the sealer it is recommended to use Isopropyl alcohol to remove remains of non-polymerized preparation.

Product setting time is about 3 hours and it complies with the standard EN ISO 6876 Root canal sealing materials. Its shortening could result in polymerization shrinkage and finally cause leaking filling.

If re-end treatment is needed removing of the sealer is possible using Canal Clean or Methylene Chloride.

CONTRAINdications

Do not use SYNTEX in patients with hypersensitivity to the product components.

PRECAUTIONS

Use a dental dam while working with the product. Use protective gloves, goggles and clothes. Causes skin irritation. Avoid skin and oral mucosa contact with non-polymerized material. Harmful if swallowed. Irritating to oral mucosa. In case of contact, immediately flush with large amounts of water. Product irritating to eyes. In case of contact, immediately flush with plenty of water. Do not induce vomiting if swallowed. Drink large amounts of water. Seek medical advice. Avoid pushing the product through the apex as it can cause irritation of the periapical tissues. It is recommended to work under X-ray and endometer guidance.

STORAGE

Store in a sealed package in a ventilated room at temperature of 2-25 °C. Protect from sunlight. Keep out of reach of children. The expiry is placed on the product label.

Once the container is open, the expiry date doesn't change if it is securely closed after each use.

When closing the syringe avoid mixing the pastes from the syringe with the pastes remaining on the closing cap (white with yellow). It can result in material polymerization and opening of the syringe will be impossible.

HANDLING THE PRODUCT CONTAINER

Empty containers should be disposed of responsibly or returned to the manufacturer.

PACKAGING

Syringe of 10 g of the preparation, mixing tips, mixing papers.

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