3M ESPE
Ketac™ Fil Plus
Glass Ionomer Filling Material
Glasionomer Füllungsmaterial
Matériau d'obturation vitro-ionomère
Materiale di riempimento a base di vetroionomero
Material de relleno de ionómero vítreo
Material de obturação de ionómero vítreo
Glasionomeer vulmateriaal
Γαλοϊονομερες εμφρακτικό υλικό
Glasjonomer-fyllningsmaterial
Lasi-ionomeeri-täyttemateriaali
Glasionomer-fylndningsmateriale
Glassionomer fyllingsmateriale

Instructions for Use
Gebrauchsinformation
Mode d'emploi
Informazioni per l'uso
Información de uso
Instruções de uso

Gebruiksformatie
Πληροφορίες χρήσεως
Bruksanvisning
Käyttöinfo
Brugsanvisning
Bruksinformasjon

	
3M ESPE AG	CE 0123
Dental Products <p>8-82229 Seefeld - Germany</p>	
3M ESPE	
Dental Products <p>St. Paul, MN U.S.A. 55144-1000</p>	
3M ESPE Customer Care/MSDS Information in U.S.A. and Canada 1-800-634-2249.	
3M, ESPE and Ketac are trademarks of 3M or 3M ESPE AG. © 3M 2008. All rights reserved.	

ENGLISH

Product Description
Ketac™ Fil Plus, manufactured by 3M ESPE, is a glass ionomer filling material for manual mixing. The following colors are available: A1/bright, A2/bright-yellow, A3/yellow, A3.5/brown, A4/yellow-brown, C3/gray-brown.

☞ For details on all products mentioned below please refer to the corresponding Instructions for Use.
Instructions for Use should not be discarded for the duration of product use.

Areas of Application

- Class III and V fillings
- Wedge-shaped defects
- Class I small fillings
- Primary tooth fillings
- Restorative fillings prior to crown preparation
- Fissure sealing

Precautionary Measures

3M ESPE MSDSs can be obtained at www.mmm.com or from your local subsidiary.

Preparation

- Remove the carious substance only, no undercuts required.
- Do not prepare thin edges since the material requires a wall thickness of 0.5 mm.

Pulp Protection

The glass ionomer cement should not be applied directly to the dentin in close proximity to the pulp or the exposed pulp. Loose abrasion material in deep cavities resulting from preparation procedures should be removed with a water spray, and dried with air. Subsequently, cover areas in close proximity to the pulp with a hard-setting calcium hydroxide material. If a sufficient dentin layer remains in the area of the pulp, no lining needs to be applied.

Conditioning

► For optimal chemical bonding to the teeth the smear layer produced in the preparatory steps must be carefully removed by applying Ketac Conditioner, manufactured by 3M ESPE, to the prepared surfaces and allowing the substance to react for 10 sec.

- Then rinse with copious amounts of water.
- Blow the cavity dry in only 2-3 short intervals with water and oil-free air or dab dry with cotton pellets. Do not over-dry!
 - The cavity should be just so dry that the surface has a matt shiny appearance.** Excessive drying can result in post-operative sensitivity after filling.
- Avoid further contamination.

Dosing

- Shake bottle to loosen the powder.
- Use a mixing ratio [weight ratio] of 3.2 parts powder (1 spoonful with leveled powder surface) : 1 part liquid (1 drop).
- To remove the powder, wipe off the spoon level at the plastic insert. Avoid compressing the powder.
- Dose powder and liquid in adjacent spots on the same mixing block.
 - Hold bottle upright when dosing the liquid.
 - The drop dosing device must be free of any dried residue of the liquid.
- The powder and liquid bottles must be carefully re-sealed after use.

Mixing

- Mix and process Ketac Fil Plus at 20–25°C/68–77°F (ambient temperature).
- Use a metal or plastic cement spatula and a mixing block or glass plate for mixing.
- Ideally, the powder should be transferred into the liquid in two aliquots.
- The paste should be streaked out repeatedly until a homogenous consistency is achieved.

Application

- Apply the Ketac Fil Plus in several increments.
- Carefully avoid contamination during application.
- The material sticks to metal instruments and should be washed off with cold water prior to setting.

Times

At an ambient temperature of 23°C/73°F and 50% relative humidity the following times apply:

	min:sec
Mixing	1:00
Processing from start of mixing	2:00
Setting from start of mixing	7:00

At higher temperatures, the time available for processing is shortened, whereas it is prolonged at lower temperatures (e. g. mixing on a chilled glass plate). The processing time is also shortened if the amount of powder exceeds the recommendation. Exceeding the indicated processing time may lead to reduced bonding to enamel and dentin.

Filling Protection

In the first phase of its 2-phase setting process, glass ionomers are sensitive to moisture. Subsequently, there is the possibility of dehydration which may be prevented by the use of Ketac™ Glaze, manufactured by 3M ESPE, directly after shaping.

- Apply Ketac Glaze to all exposed surfaces of the glass ionomer cement and cure with light for 10 sec.

Finishing

- Beginning no earlier than 7 minutes after the start of mixing, use Arkansas stone, fine-grain diamonds, polishing discs of decreasing grain size or silicone polishers for finishing.
- If desired, apply Ketac Glaze again for surface sealing and cure with light for 10 sec.

Storage and Shelf-life

Store the product at 15–25°C/59–77°F. Do not use after the expiry date.

Customer Information

No person is authorized to provide any information which deviates from the information provided in this instruction sheet.

Warranty

3M ESPE warrants this product will be free from defects in material and manufacture. 3M ESPE MAKES NO OTHER WARRANTIES INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and 3M ESPE's sole obligation shall be repair or replacement of the 3M ESPE product.

Limitation of Liability

Except where prohibited by law, 3M ESPE will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

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