**Product : ProTaper® Universal (Treatment)****Ref doc. : F19 02 85.X / 01 / 2006 updated 05/2006****CE  
0459**

		Page
GB	DIRECTION FOR USE	2
F	PROTOCOLE D'UTILISATION	6
D	GEBRAUCHSANWEISUNG	10
E	ISTRUCCIONES DE USO	14
I	ISTRUZIONI PER L'USO	18
DK	BRUGSANVISNING	22
S	BRUKSANVISNING	26
FIN	KÄYTTÖOHJEET	29
NL	GEbruiksvoorschrift	33
P	PROTOCOLO DE UTILIZAÇÃO	37
GR	ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ ΓΙΑ	41

**EC REP**Dentsply France SAS  
4, rue Michael Faraday  
F – 78180 Montigny Le Bretonneux

## Directions for use ProTaper® Universal

### ProTaper® instruments for endodontic treatment:

- ProTaper® shaping files
- ProTaper® finishing files

**COMPOSITION:** the cutting part of these instruments is made of a nickel-titanium alloy.

### 1. INDICATIONS FOR USE:

These instruments are to be used only in a clinical or hospital environment, by qualified users.

**Application field:** shaping and cleaning of the root canal system.

### 2. CONTRAINDICATIONS:

None known.

### 3. WARNINGS:

None known.

### 4. PRECAUTIONS:

- Multiple use disinfection and resterilization cycles may lead to increased risk of file separation.
- These instruments should not be immersed in a sodium hypochlorite solution.
- Instrument decontamination: strictly follow decontamination instructions from the manufacturer.
- Irrigate abundantly and frequently
- Establish a reproducible glide path using hand files, at least an Iso 015 size.
- Use in a constant rotation at a speed of 150-350 rpm with light apical pressure.
- Clean flutes frequently and check for signs of distortion or wear.
- For optimal usage, torque control devices are recommended.
- Use the Shaping Files (S1, S2 and SX) with a brushing action on the withdrawal stroke in order to create straight line radicular access.
- Use the Finishing Files (F1, F2, F3, F4 and F5) with no brushing action.
- Use the appropriate finishing files to passively follow the canal to the working length, and then withdraw immediately.

### 5. ADVERSE REACTIONS:

In the present technical state, no adverse reaction has been reported so far.

### 6. STEP BY STEP INSTRUCTIONS FOR PROTAPER® FILES:

- 1) Create straight line access to canal orifice
- 2) Always irrigate and confirm a reproducible glide path with an ISO 015 hand file.
- 3) Protocol of use:
  - Locate the orifice.
  - Use passively an ISO 015 hand file to resistance.
  - Use Shaping File S1 with brushing action until the same depth as the ISO 015 hand file is reached.
  - Repeat this sequence until the working length is determined with an ISO 015 hand file and reached with S1.
  - Use shaping file S2 with brushing action until the working length is reached.
  - reconfirm working length.
  - Use finishing file F1 (non-brushing motion) with each insertion deeper than the previous insertion until the working length is reached.
  - Gauge the foramen with hand files.
  - Use the appropriate finishing file (F2, F3, F4, F5) with the same non-brushing motion to working length if additional enlargement is required or if foramen is larger.

If necessary, use the SX with a brushing motion to move the coronal aspect of the canal away from the furcation and / or to create more coronal shape.

## **7. DISINFECTION, CLEANING AND STERILIZATION:**

First clean and decontaminate the instruments and their supports. Then decontaminate and sterilize them out of their environment (in bags) with the autoclave at 134°C and 3 bar during 18 minutes.

### **7.1 General recommendation**

- Single use marked instruments are not approved for re-use.
- The user is responsible for the sterility of the product for the first cycle and each further usage as well as for the usage of damaged or dirty instruments.
- For your own safety, please wear personal protective equipments (gloves, glasses).
- Use only disinfecting solution which is approved for their efficacy (DGHM-listing, CE marking, FDAapproval)
- Tungsten Carbide burs, Plastic Support, Hand instruments and NiTi instruments are degraded by Peroxide hydrogen (H<sub>2</sub>O<sub>2</sub>) solution.
- NiTi Instruments are degraded if immersed more than 5 minutes in a solution of NaOCl at more than 5 %.
- Instruments made of aluminium are degraded in presence of caustic soda solutions with mercury salt.
- Do not use acid (pH < 6) or alkaline (pH > 8) solutions.
- After 5 cycles, degradation of marking may appear.

## 7.2 Step-by-step procedure

### A. Devices

			A2. Implantable devices			
			A1. Instruments	Following usages		
				First usage		
	Operation	Operating mode	Warning			
1.	Pre-Disinfection or Decontamination	Soak immediately just after usage all instruments in a disinfectant solution combined with proteolytic enzyme if possible.	<ul style="list-style-type: none"> <li>- Follow instructions and observe concentrations and immersion time given by the manufacturer (an excessive concentration may cause corrosion or others defects on instruments).</li> <li>- The disinfectant solution should be aldehyde free (to avoid blood impurities fixation).</li> <li>- Do not use disinfectant solution containing Phenol or any products which are not compatible with the instruments (See general recommendation).</li> <li>- For visible impurities are observed on instruments a pre-cleaning is recommended by brushing them manually with soft material.</li> </ul>		X	
2a.	Automated Cleaning / Disinfection	<ul style="list-style-type: none"> <li>- Disassemble the devices (silicone stops to be removed)</li> <li>- Place them in a kit, support or container</li> <li>- Put them in the thermodisinfectant (for at least 10 min at 93°C or Ao value &gt; 3000).</li> </ul>	<ul style="list-style-type: none"> <li>- Discard any instruments with large obvious defects (broken, bent).</li> <li>- Avoid any contact between instruments or posts when placing in the thermodisinfectant.</li> <li>- Follow instructions and observe concentrations given by the manufacturer (see also general recommendation).</li> <li>- The thermodisinfectant is not recommended for instruments made of aluminium, tungsten carbide or carbon steel. In case of chemical disinfection, there is a danger of remnants of the disinfectant on the instruments.</li> <li>- Sufficient rinsing step should be available in the program. Use purified water for rinsing (max 10 germs/ml and max 0.25 endotoxin units/ml) like Aqua Purificata.</li> <li>- Use filtered air for drying the instruments or posts.</li> <li>- Use only approved thermodisinfectant according to EN ISO 15883, maintain and calibrate it regularly.</li> <li>- If possible, an automated procedure is preferred.</li> </ul>	X	X	X
2b.	Manual Cleaning / Disinfection	<ul style="list-style-type: none"> <li>- Disassemble the devices (silicone stops to be removed)</li> <li>- Immerse them in the disinfectant solution assisted with ultrasonic device if suitable.</li> <li>- Rinse the devices thoroughly with clean, demineralised or distilled water and then dry them with filtered compressed air.</li> </ul>	<ul style="list-style-type: none"> <li>- No visible impurities should be observed on the instruments.</li> <li>- Discard any instruments with large obvious defects (broken, bent, and twisted).</li> <li>- Avoid any contact between instruments or posts when placing in the solution, use kits, supports or container.</li> <li>- Follow instructions and observe concentrations and time given by the manufacturer (see also general recommendation).</li> <li>- Use purified water for rinsing (max 10 germs/ml and max 0.25 endotoxin units/ml) like Aqua Purificata.</li> <li>- If a disinfectant solution contains a corrosion inhibitor, it is recommended to rinse the instruments just before autoclaving.</li> </ul>	X	X	X
3.	Inspection	<ul style="list-style-type: none"> <li>- Inspect devices and sort out those with defects.</li> <li>- Assemble the devices (stops)</li> </ul>	<ul style="list-style-type: none"> <li>- Dirty instruments must be cleaned and disinfected again.</li> <li>- Discard instruments which show any deformations (bent, twisted), damages (broken, corroded) or defects (loss of colour coding or marking) affecting the resistance, the safety or the performance of the instrument or posts.</li> <li>- Protect carbon steel bur with corrosion inhibitor before packaging.</li> </ul>		X	X
4.	Packaging	- Pack the devices in "Sterilisation pouches".	<ul style="list-style-type: none"> <li>- Check the validity period of the pouch given by the manufacturer to determine the shelf life.</li> <li>- Use packaging which are resistant up to a temperature of 141°C (286°F) and in accordance with EN ISO 11607.</li> </ul>	X	X	X
5.	Sterilization	<ul style="list-style-type: none"> <li>- Steam sterilisation at: 134°C (237°F) during 18 min for endodontics instruments and root canal filling instruments</li> <li>- 134 °C (237°F) during at least 5 min for all other devices.</li> </ul>	<ul style="list-style-type: none"> <li>- The instruments, posts and the plastic supports must be sterilized according to the packaging labelling.</li> <li>- Use fractionated vacuum or gravity (less preferred) autoclaves (according to EN 13060, EN 285).</li> <li>- Use validated sterilisation procedure according to ISO 17665-1</li> <li>- Respect maintenance procedure of the autoclave device given by the manufacturer.</li> <li>- Use only the listed sterilization procedures.</li> </ul>	X	X	X
6.	Storage	- Keep devices in sterilization packaging in a dry and clean environment.	- Sterility cannot be guaranteed if packaging is open, damaged or wet (check the packaging before using the instruments).	X	X	X

## B. Filling material

	Operation	Operating mode	Warning
1.	Disinfection	- Immerse the obturation devices in alcohol during 30 sec at ambient temperature.	- Do not use disinfectant solution containing Phenol or any products which are not compatible with the instruments (See general recommendation).

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