General Information:

Adper™ Scotchbond™ 1 XT Adhesive, manufactured by 3M ESPE, is a simple, moist bonding adhesive containing 10%, 5nm colloidal filler.

Adper Scotchbond 1 XT adhesive offers the dental practitioner a wide range of applications. These include bonding to all classes of direct composite restorations as well as procedures involving porcelain, composite, metal repair, set amalgam, root surface desensitization and bonding of porcelain veneers with RelyX™ Venerer Cement System and RelyX™ Ceramic Primer, manufactured by 3M ESPE.

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Precautions for Dental Personnel and Patients

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Adper Scotchbond 1 XT adhesive contains acrylates including HEMA (2-hydroxyethylmethacrylate). Avoid use of this product on patients with known acrylate allergies. To reduce the risk of allergic response, minimize exposure to these materials. In particular, avoid exposure to uncured resins. Use of protective gloves and a sealach technique is recommended. If skin contact occurs, wash skin with soap and water. Acrylates may penetrate commonly-used gloves. If adhesive contacts glove, remove and discard glove, wash hands immediately with soap and water and then re-glove. If accidental contact with eyes or prolonged contact with oral soft tissue occurs, flush with large amounts of water. If irritation persists, consult a physician. See Material Safety Data Sheet (MSDS) for additional cautionary information. You may obtain the current MSDS by visiting the website: <http://www.3m.com/MSDS> or contacting your 3M ESPE Dental Products representative.

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Etchant Syringe Assembly:

- Protective eyewear for patients and staff is recommended when using the delivery system.
- Prepare the delivery system: Remove cap from etchant syringe and SAVE. Twist a large disposable tip securely onto the syringe. Holding the syringe with the tip away from the patient and any dental staff, express a small amount of etchant onto a dispensing pad or a 2 x 2 pad to assure that the delivery system is not clogged. If clogged, remove the dispensing tip and express a small amount of etchant directly from the syringe. Remove any visible plug, if present, from the syringe opening. Replace dispensing tip and again express etchant. If clog remains, discard dispensing tip and replace with a new one. Bend the dispensing tip to a desired angle. Place bend midway along tip. Do not bend dispensing tip at its hub as this may cause the tip to break free.
- Delivery syringe setup: Remove used dispensing tip and discard. Twist on storage cap. **Storage of the delivery syringe with a used dispensing tip or without storage cap will allow drying of the etchant and consequent clogging of the system.**
- Replace storage cap with a new dispensing tip at next use.
- If desired, the etchant may be extruded onto a dispensing pad and applied with a brush or other appropriate instrument.
- If a liquid etchant is desired, the etchant may be dispensed into a dappen dish and stirred to increase its fluidity.

- Discinfection: Discard used dispensing tip. Replace syringe cap. Discinfect the capped syringe in the same manner as nonmiserable handicpes, air/water syringe and ultrasonic scales following American Dental Association (ADA) and Centers for Disease Control (CDC) recommendations. (Council on Dental Materials, Instruments and Equipment and Council on Dental Therapeutics. Infection control recommendations for the dental office and the dental laboratory. JADA 116(2):241-248, 1988.)

Dispensing Adper Scotchbond 1 XT Adhesive

Unit Dose:

Attention: To minimize risk of accidental contact with eyes and skin, hold the shaft of the disposable applicator over the opening where the applicator enters the tip package with your thumb and index finger. Do not activate the unit dose without a disposable applicator. With your thumb and index finger of the other hand, squeeze the large blister to transfer the adhesive into the chamber enclosing the applicator. Briefly spin the applicator to fully saturate with adhesive.

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Instructions for direct light cure restorations in enamel and dentin:

- Isolation: Rubber dam is the preferred method of isolation.
- Cavity preparation: Prepare cavity with minimal tooth reduction. Bevel cavosurface enamel margins.
- Etching: Apply Scotchbond etchant to enamel and dentin. Wait 15 seconds. Rinse for 10 seconds. Blot excess water using a cotton pellet or mini-sponge. The surface should appear glistening without pooling of water.
- Adhesive: Immediately after blotting, apply 2-3 consecutive coats of adhesive to etched enamel and dentin for 15 seconds with gentle agitation using a fully saturated applicator. Gently air thin for five seconds to evaporate solvents. Light cure for 10 seconds.
- Restorative placement, cure and finishing: Refer to manufacturer's instructions for placement, cure and finishing of restorative materials.

Instructions for bonding porcelain veneers:

- Silane treatment: Porcelain bonding surfaces should have been etched using hydrofluoric acid by the dental laboratory. Apply RelyX Ceramic Primer (No. 2721) to the bonding surface of the veneer. Dry for 5 seconds.
- Clean the prepared teeth in preparation for seating and bonding using a plain flour of pumice slurry. Rinse thoroughly and dry.
- Try in veneer with RelyX™ Try-In Paste, manufactured by 3M ESPE. After try in, isolate from adjacent teeth with wax matrix strip.
- Etching: Apply Scotchbond etchant to both enamel and dentin. Wait 15 seconds. Rinse for 10 seconds. Blot excess water using a cotton pellet or mini-sponge. The surface should appear glistening without pooling of water. An additional 15 second etch time may be appropriate for teeth that were not prepared using a diamond or bur. Residual organic matter can also require additional etch time.
- Adhesive: Immediately after blotting, apply 2-3 consecutive coats of adhesive to etched enamel and dentin for 15 seconds with gentle agitation using a fully saturated applicator. Gently air thin for five seconds to evaporate solvents. Do not light cure.
- Adhesive application to veneer: Apply 1 coat of adhesive to the acid etched, silane treated veneer. Dry thoroughly. Do not light cure.
- Luting material application to veneer: Apply RelyX Venerer Cement to the bonding surface of veneer.

- Seating and curing: Carefully seat the veneer. Clean excess luting cement from the veneer margins. Cure each area of the veneer for times recommended by the luting cement manufacturer. We recommend curing the initial margin first, followed by the body and the incisal margin. Avoid direct contact with the light-gauge during initial curing.

Instructions for bonding to composite and set amalgam:

- Isolation: Rubber dam is the preferred method of isolation.
- Roughen the existing material: Roughen the surface of existing amalgam or composite using either a bur, diamond or a sandblast technique.
- Etching: Apply Scotchbond etchant to enamel, dentin and existing restorative material. Wait 15 seconds. Rinse for 10 seconds. Blot excess water using a cotton pellet or mini-sponge. The surface should appear glistening without pooling of water.

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5) Masking: In the case of set amalgam, mask the metal surface with a thin layer (0.25 – 0.5mm) of the appropriate 3M ESPE Masking Agent shade, manufactured by 3M ESPE, with a brush. Light cure each layer for 20 seconds.

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Instructions for porcelain repair:

- Isolation: Rubber dam is the preferred method of isolation.
- Preparation: Clean the surface to be repaired with a slurry of plain flour of pumice. Rinse and dry thoroughly. Roughen the surface of existing metal or porcelain using either a bur, diamond or a sandblast technique. Be careful to remove all loose porcelain and bevel the margin. Remove surface glaze 1mm beyond the margin.
- Etching: Apply Scotchbond etchant to all substrates. Wait for 15 seconds. Rinse for 10 seconds. Dry 5 seconds.
- Silane treatment for porcelain and metal: Apply RelyX Ceramic Primer (No. 2721) to the etched surface and dry.
- Adhesive: Apply 2 consecutive coats of adhesive to etched enamel 1 XT adhesive to silane treated porcelain or metal. Dry gently for 5 seconds. Light cure for 10 seconds.
- Masking: To opacity the metal before the final composite placement, mask the metal surface with a thin layer (0.25 – 0.5mm) of the appropriate 3M ESPE Masking Agent shade using a brush. Light cure each layer for 20 seconds.

- Restorative placement, cure and finishing: Refer to manufacturer's instructions for placement, cure and finishing of restorative material.

Instructions for root surface desensitization:

- Lightly clean the root surface with flour of pumice. Rinse and blot dry.
- Etching: Apply Scotchbond etchant to the etched dentin. Wait for 15 seconds. Rinse for 10 seconds. Blot excess water using a cotton pellet or mini-sponge. The surface should appear glistening without pooling of water.
- Adhesive: Immediately after blotting, apply 2-3 consecutive coats of adhesive to etched enamel and dentin for 15 seconds with gentle agitation using a fully saturated applicator. Gently air thin for five seconds to evaporate solvents. Light cure for 10 seconds. Apply 2 additional coats of adhesive. Dry gently for 5 seconds. Light cure for 10 seconds.
- Remove the oxygen inhibited layer with a moistened gauze.

Instructions for Amalgam and Indirect Bonding Procedures (using RelyX ARC Composite Resin cement):

Physical properties of today's esthetic indirect restorations require that they be bonded into place to maximize the strength of the restoration as well as the tooth. A general perception may exist that light cured adhesives cannot be used for indirect restorations. It's true that many conventional light cured adhesives have a higher film thickness and can not be used under a fixed prosthesis. However, Adper Scotchbond 1 XT adhesive is ethanol/water based, has a low film thickness (approximately 10µm) and should not interfere with the seating of indirect restorations.

Note:

Care is required with any bonding agent used beneath precision castings because added film thickness may preclude accurate seating. Avoid adhesive pooling in areas of the preparation that would offset the fit of any prosthetic device.

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- Remove temporary restoration. Trial-fit the final restoration with light finger pressure to evaluate the fit, shade and marginal integrity. Adjust if necessary.
- Prepare the bonding surface of the indirect restoration and the core build up, if applicable. Porcelain bonding surfaces should have been etched with hydrofluoric acid by the dental laboratory. Metal and amalgam bonding surfaces should be roughened, preferably using an air abrasion system, diamond or bur. Any composite surfaces should be roughened with a diamond, bur or air abrasion system. Glass ionomer build-ups should be pumiced with a slurry of plain flour of pumice.
- Silane treatment (porcelain or porcelain/metal indirect restorations): Apply RelyX Ceramic Primer to the bonding surface of the indirect restoration. Dry for 5 seconds.
- Clean the prepared teeth in preparation for seating and bonding using a plain flour of pumice slurry. Rinse and dry thoroughly, isolate from moisture and adjacent teeth.
- Etching: Apply Scotchbond etchant to both enamel and dentin. Wait 15 seconds. Rinse for 10 seconds. Blot excess water using a cotton pellet or mini-sponge. The surface should appear glistening without pooling of water.
- Adhesive: Immediately after blotting, apply 2-3 consecutive coats of adhesive to etched enamel and dentin for 15 seconds with gentle agitation using a fully saturated applicator. Gently air thin for five seconds to evaporate solvents being careful to avoid excess adhesive on all prepared surfaces.
- Light cure preparation for 10 seconds per bonding surface.
- Dispense the appropriate amount of cement onto a mixing pad and mix for 10 seconds.
- Apply and evenly distribute a thin layer of cement to the bonding surface of the indirect restoration.
- Slowly seal and hold restoration in proper occlusion. Begin clean-up of excess cement approximately 3-5 minutes after seating. **Optional: If excess cement is removed immediately after seating, each cement surface/margin must be light cured for 40 seconds.**

- Once the restoration is seated, each cement surface/margin may be light cured for 40 seconds or allowed to self cure for 10 minutes. **Note: for porcelain and pre-cured composite restorations, each cement surface/margin must be light cured for 40 seconds.**
- Instruct patient to avoid applying any pressure for 10-15 minutes.

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- Prepare the endodontically treated tooth to receive the post (a root apex sealer and gutta percha filling approximately one third of the root canal are recommended). **Trial fit and adjust post as needed.** Bond to cast posts can be enhanced by using an air abrasion system and then applying RelyX Ceramic Primer. Dry for 5 seconds.
- Etching: Apply Scotchbond etchant to the prepared tooth. Wait 15 seconds. Rinse for 10 seconds. Dry for 2 seconds. Remove excess moisture with an absorbent paper point.
- Adhesive: Apply a uniform coat to etched enamel and dentin. Remove excess pooled adhesive with absorbent paper point. Air thin for five seconds to evaporate solvents.
- Light cure for 10 seconds. (A light transmitting post of appropriate size may be used for additional curing.)
- Dispense the appropriate amount of cement onto a mixing pad and mix for 10 seconds.
- Apply cement to the bonding surface of the preparation (in and around canal using a periodontal probe). Place a thin layer of mixed cement on post.
- Seat the post. While holding in place remove excess cement. Light cure for 40 seconds from the occlusal surface to allow immediate placement of core build-up material

Instructions for bonding amalgam to both structure:

- Isolation: Rubber dam is the recommended method of isolation.
- Cavity preparation: Prepare a standard amalgam cavity preparation. Roughen residual restorative materials with an air abrasion system or a bur.
- Matrix application: Lightly lubricate the inner surface of the matrix band with hard wax or petroleum jelly before placement.
- Etching: Apply Scotchbond etchant to enamel, dentin and any residual restorative. Wait for 15 seconds. Rinse for 10 seconds. Blot excess water using a cotton pellet or mini-sponge. The surface should appear glistening without pooling of water.
- Adhesive: Immediately after blotting, apply 2-3 consecutive coats of adhesive to etched enamel, dentin and any residual restorative material for 15 seconds with gentle agitation using a fully saturated applicator. Gently air thin for five seconds to evaporate solvents.
- Light cure for 10 seconds.
- Dispense the appropriate amount of cement onto a mixing pad and mix for 10 seconds.
- Use a brush or appropriate applicator to place cement in adhesive-sealed preparation.

Trisilane amalgam preparation of cement.

- Condense and burnish amalgam in the usual way.

- Instruct patient to avoid applying any pressure for 10-15 minutes.

Additional notes:

- Brush handles can be disinfected in the same manner as nonmiserable handicpes, air/water syringe and ultrasonic scales following American Dental Association (ADA) and Center for Disease Control (CDC) recommendations. (Council on Dental Materials, Instruments and Equipment and Council on Dental Therapeutics. Infection control recommendations for the dental office and the dental laboratory. JADA 116(2):241-248, 1988.)

Storage and use:

- Adper Scotchbond 1 XT system can be stored at room temperature.
- Adper Scotchbond 1 XT adhesive should be capped immediately after use to minimize evaporation.
- Do not expose materials to elevated temperature or intense light.
- Do not store products in proximity to eugenol containing products.
- This system is designed to be used at room temperature of approximately 21-24°C or 70-75°F.
- Shelf life of the unit dose is 24 months at room temperature. Shelf life of the vial delivery is 36 months at room temperature. See outer package for expiry date.
- For cleanup, Scotchbond etchant can be removed with water, while the uncured adhesive can be removed with alcohol. No person is authorized to provide any information which deviates from the information provided in this instruction sheet.

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## ENGLISH

### General Information:

Adper™ Scotchbond™ 1 XT Adhesive, manufactured by 3M ESPE, is a simple, moist bonding adhesive containing 10%, 5mm colloidal filler.

Adper Scotchbond 1 XT adhesive offers the dental practitioner a wide range of applications. These include bonding to all classes of direct composite restorations as well as procedures involving porcelain, composite, metal repair, set amalgam, root surface desensitization and bonding of porcelain veneers with RelyX™ Veneer Cement System and RelyX™ Ceramic Primer, manufactured by 3M ESPE.

After light curing Adper Scotchbond 1 XT adhesive, it may also be used for amalgam and indirect bonding procedures when combined with RelyX™ ARC adhesive resin cement, manufactured by 3M ESPE. Compatibility with indirect bonding procedures is due to the low film thickness (approximately 10µm) of cured Adper Scotchbond 1 XT.

Adper Scotchbond 1 XT adhesive is available in two delivery systems, a unit dose delivery and a multi-use vial dispenser.

Use of etchant is critical for both enamel and dentin surfaces.

### Recommendations:

Use Vitrebond™ Light Cure Glass Ionomer Liner/Base, manufactured by 3M ESPE, in areas of deep cavity excavation such as Class I and II restorations. If pulp exposure has occurred, use a minimum amount of calcium hydroxide followed by an application of Vitrebond liner/base. Adper Scotchbond 1 XT adhesive will bond to Vitrebond liner/base whether or not the ionomer was treated with etchant.

Adper Scotchbond 1 XT adhesive includes an etch of enamel and dentin as a part of the procedure. It is recommended that the surfaces be left moist after rinsing. Excess surface moisture should be removed by blotting.

Adper Scotchbond 1 XT adhesive is cured by exposure to visible light. The light curing times instructed with this product assumes the use of a 3M ESPE light curing unit, manufactured by 3M ESPE, or other dental visible curing light of comparable intensity. Curing lights should be checked often for proper output using a reliable light metering system.

Air used for drying should be free of oil and water contaminants.

### Precautions for Dental Personnel and Patients

Scotchbond™ Etchant, manufactured by 3M ESPE, contains 35 weight % phosphoric acid.

Protective eyewear for patients and dental staff is recommended when using etchants. Avoid contact with oral soft tissue, eyes and skin. If accidental eye contact occurs, flush immediately with large amounts of water. For eye contact consult a physician.

Adper Scotchbond 1 XT adhesive contains acrylates including HEMA (2-hydroxyethylmethacrylate). Avoid use of this product on patients with known acrylate allergies. To reduce the risk of allergic response, minimize exposure to these materials. In particular, avoid exposure to uncured resins. Use of protective gloves and a no-touch technique is recommended. If skin contact occurs, wash skin with soap and water. Acrylates may penetrate commonly-used gloves. If adhesive contacts glove, remove and discard glove, wash hands immediately with soap and water and then re-glove. If accidental contact with eyes or prolonged contact with oral soft tissue occurs, flush with large amounts of water. If irritation persists, consult a physician. See Material Safety Data Sheet (MSDS) for additional cautionary information. You may obtain the current MSDS by visiting the website: <http://www.3m.com/MSDS> or contacting your 3M ESPE Dental Products representative.

### Sensitivity:

Some patients may experience transitory postoperative sensitivity. The risk of sensitivity can be minimized by the following measures:

Tooth Preparation:

Remove minimal tooth structure.

Use proper isolation. Use of a rubber dam is highly recommended. Use adequate pulp protection. Use a glass ionomer or resin-modified glass ionomer liner/base (3M ESPE Vitrebond) in areas of deep excavation.

Adhesive Application:

Use of compressed air is not recommended to remove pooled water remaining after the etch step–blot excess moisture from the preparation using a cotton pellet or mini-sponge.

Apply adhesive immediately after blotting.

Restorative:

Place restorative material in increments, curing each increment separately.

Adequately cure restorative according to instructions for shade and thickness of restorative and light exposure time.

Adjust occlusion carefully. Check for hyperocclusion, particularly in lateral excursion contacts.

### Etchant Syringe Assembly:

1) Protective eyewear for patients and staff is recommended when using the delivery system.

2) Prepare the delivery system: Remove cap from etchant syringe and SAVE. Twist a blue disposable tip securely onto the syringe. Holding the syringe with the tip away from the patient and any dental staff, express a small amount of etchant onto a dispensing pad or a 2 x 2 gauze to assure that the delivery system is not clogged. If clogged, remove the dispensing tip and express a small amount of etchant directly from the syringe. Remove any visible plug, if present, from the syringe opening. Replace dispensing tip and again express etchant. If clog remains, discard dispensing tip and replace with a new one. Bend the dispensing tip to a desired angle. Place bend midway along tip. Do not bend dispensing tip at its hub as this may cause the tip to break free.

3) Delivery system storage: Remove used dispensing tip and discard. Twist on storage cap. Storage of the delivery syringe with a used dispensing tip or without storage cap will allow drying of the etchant and consequent clogging of the system. Replace storage cap with a new dispensing tip at next use.

4) If desired, the etchant may be extruded onto a dispensing pad and applied with a brush or other appropriate instrument.

5) If a liquid etchant is desired, the etchant may be dispensed into a dappen dish and stirred to increase its fluidity.

6) Disinfection: Discard used dispensing tip. Replace syringe cap. Disinfect the capped syringe in the same manner as nonimmersible handpieces, air/water syringe and ultrasonic scalers following American Dental Association (ADA) and Centers for Disease Control (CDC) recommendations. (Council on Dental Materials, Instruments and Equipment and Council on Dental Therapeutics. Infection control recommendations for the dental office and the dental laboratory. JADA 116(2):241-248, 1988.)

### Dispensing Adper Scotchbond 1 XT Adhesive

Unit Dose:

Attention: To minimize risk of accidental contact with eyes and skin, hold the shaft of the disposable applicator over the opening where the applicator enters the foil package with your thumb and index finger. Do not activate the unit dose without a disposable applicator. With your thumb and index finger of the other hand, squeeze the large blister to transfer the adhesive into the chamber enclosing the applicator. Briefly spin the applicator to fully saturate with adhesive.

Vial:

Pinch the sides of the cap to release the locking mechanism and flip the cap back to reveal the dispensing tip. Squeeze out the exact number of drops you need into the disposable mixing well. When finished, flip the cap back until it is secured by its locking mechanism.

### Instructions for direct light cure restorations in enamel and dentin:

1) Isolation: Rubber dam is the preferred method of isolation.

2) Cavity preparation: Prepare cavity with minimal tooth reduction. Bevel cavosurface enamel margins.

3) Etching: Apply Scotchbond etchant to enamel and dentin. Wait 15 seconds. Rinse for 10 seconds. Blot excess water using a cotton pellet or mini-sponge. The surface should appear glistening without pooling of water.

4) Adhesive: Immediately after blotting, apply 2-3 consecutive coats of adhesive to etched enamel and dentin for 15 seconds with gentle agitation using a fully saturated applicator. Gently air thin for five seconds to evaporate solvents. Light cure for 10 seconds.

5) Restorative placement, cure and finishing: Refer to manufacturer's instructions for placement, cure and finishing of restorative materials.

### Instructions for bonding porcelain veneers:

1) Silane treatment: Porcelain bonding surfaces should have been etched using hydrofluoric acid by the dental laboratory. Apply RelyX Ceramic Primer (No. 2721) to the bonding surface of the veneer.

Dry for 5 seconds.

2) Clean the prepared teeth in preparation for seating and bonding using a plain flour of pumice slurry. Rinse thoroughly and dry.

3) Try in veneer with RelyX™ Try-In Paste, manufactured by 3M ESPE. After try in, isolate from adjacent teeth with clear matrix strip.

4) Etching: Apply Scotchbond etchant to both enamel and dentin. Wait 15 seconds. Rinse for 10 seconds. Blot excess water using a cotton pellet or mini-sponge. The surface should appear glistening without pooling of water. An additional 15 second etch time may be appropriate for teeth that were not prepared using a diamond or bur. Residual organic matter can also require additional etch time.

5) Adhesive: Immediately after blotting, apply 2-3 consecutive coats of adhesive to etched enamel and dentin for 15 seconds with gentle agitation using a fully saturated applicator. Gently air thin for five seconds to evaporate solvents. Do not light cure.

6) Adhesive application to veneer: Apply 1 coat of adhesive to the acid etched, silane treated veneer. Dry thoroughly. Do not light cure.

7) Luting material application to veneer: Apply RelyX Veneer Cement to the bonding surface of veneer.

8) Seating and curing: Carefully seat the veneer. Clean excess luting cement from the veneer margins. Cure each area of the veneer for times recommended by the luting cement manufacturer. We recommend curing the gingival margin first, followed by the body and the incisal margin. Avoid direct contact with the light-guide during initial curing.

### Instructions for bonding to composite and set amalgam:

1) Isolation: Rubber dam is the preferred method of isolation.

2) Roughen the existing material: Roughen the surface of existing amalgam or composite using either a bur, diamond or a sandblast technique.

3) Etching: Apply Scotchbond etchant to enamel, dentin and existing restorative material. Wait 15 seconds. Rinse for 10 seconds. Blot excess water using a cotton pellet or mini-sponge. The surface should appear glistening without pooling of water.

4) Adhesive: Immediately after blotting, apply 2-3 consecutive coats of adhesive to etched enamel, dentin and existing restorative material for 15 seconds with gentle agitation using a fully saturated applicator. Gently air thin for five seconds to evaporate solvents. Light cure for 10 seconds.

5) Masking: In the case of set amalgam, mask the metal surface with a thin layer (0.25 - 0.5mm) of the appropriate 3M ESPE Masking Agent shade, manufactured by 3M ESPE, with a brush. Light cure each layer for 20 seconds.

6) Restorative placement: Refer to manufacturer's instructions for placement, cure and finishing of restorative material.

### Instructions for porcelain repair:

1) Isolation: Rubber dam is the preferred method of isolation.

2) Preparation: Clean the surface to be repaired with a slurry of plain flour of pumice. Rinse and dry thoroughly. Roughen the surface of existing metal or porcelain using either a bur, diamond or a sandblast technique. Be careful to remove all loose porcelain and bevel the margin. Remove surface glaze 1mm beyond the margin.

3) Etching: Apply Scotchbond etchant to all substrates. Wait for 15 seconds. Rinse for 10 seconds. Dry 5 seconds.

4) Silane treatment for porcelain and metal: Apply RelyX Ceramic Primer (No. 2721) to the etched surface and dry.

5) Adhesive: Apply 2 consecutive coats of Adper Scotchbond 1 XT adhesive to silane treated porcelain or metal. Dry gently for 5 seconds. Light cure for 10 seconds.

6) Masking: To specify the metal before the final composite placement, mask the metal surface with a thin layer (0.25 - 0.5mm) of the appropriate 3M ESPE Masking Agent shade using a brush. Light cure each layer for 20 seconds.

7) Restorative placement, cure and finishing: Refer to manufacturer's instructions for placement, cure and finishing of restorative material.

### Instructions for root surface desensitization:

1) Lightly clean the root surface with flour of pumice. Rinse and blot dry.

2) Etching: Apply Scotchbond etchant to the etched dentin. Wait for 15 seconds. Rinse for 10 seconds. Blot excess water using a cotton pellet or mini-sponge. The surface should appear glistening without pooling of water.

3) Adhesive: Immediately after blotting, apply 2-3 consecutive coats of adhesive to etched enamel and dentin for 15 seconds with gentle agitation using a fully saturated applicator. Gently air thin for five seconds to evaporate solvents. Light cure for 10 seconds. Apply 2 additional coats of adhesive. Dry gently for 5 seconds. Light cure for 10 seconds.

4) Remove the oxygen inhibited layer with a moistened gauze.

### Instructions for Amalgam and Indirect Bonding Procedures (using RelyX ARC adhesive resin cement):

Physical properties of today's esthetic indirect restorations require that they be bonded into place to maximize the strength of the restoration as well as the tooth. A general perception may exist that light cured adhesives cannot be used for indirect restorations. It's true that many conventional light cure adhesives have a higher film thickness and can not be used under a fixed prosthesis. However, Adper Scotchbond 1 XT adhesive is ethanol/water based, has a low film thickness (approximately 10µm) and should not interfere with the seating of indirect restorations.

### Note:

Care is required with any bonding agent used beneath precision castings because added film thickness may preclude accurate seating. Avoid adhesive pooling in areas of the preparation that would effect the fit of any prosthetic device.

### Instructions for bonding crowns, bridges (including resin-bonded bridges), inlays and onlays:

1) Remove temporary restoration. Trial-fit the final restoration with light finger pressure to evaluate the fit, shade and marginal integrity. Adjust if necessary.

2) Prepare the bonding surface of the indirect restoration and the core build up, if applicable.

Porcelain bonding surfaces should have been etched with hydrofluoric acid by the dental laboratory. Metal and amalgam bonding surfaces should be roughened, preferably using an air abrasion system, diamond or bur. Any composite surfaces should be roughened with a diamond, bur or air abrasion system. Glass ionomer build-ups should be primed with a slurry of plain flour of pumice.

3) Silane treatment (porcelain or porcelain/metal indirect restorations): Apply RelyX Ceramic Primer to the bonding surface of the indirect restoration. Dry for 5 seconds.

4) Clean the prepared teeth in preparation for seating and bonding using a plain flour of pumice slurry. Rinse and dry thoroughly, isolate from moisture and adjacent teeth.

5) Etching: Apply Scotchbond etchant to both enamel and dentin. Wait 15 seconds. Rinse for 10 seconds. Blot excess water using a cotton pellet or mini-sponge. The surface should appear glistening without pooling of water.

6) Adhesive: Immediately after blotting, apply 2-3 consecutive coats of adhesive to etched enamel and dentin for 15 seconds with gentle agitation using a fully saturated applicator. Gently air thin for five seconds to evaporate solvents being careful to avoid excess adhesive on all prepared surfaces.

7) Light cure preparation for 10 seconds per bonding surface.

8) Dispense the appropriate amount of cement onto a mixing pad and mix for 10 seconds.

9) Apply and evenly distribute a thin layer of cement to the bonding surface of the indirect restoration.

10) Slowly seat and hold restoration in proper occlusion. Begin clean-up of excess cement approximately 3-5 minutes after seating. Optional: If excess cement is removed immediately after seating, each cement surface/margin must be light cured for 40 seconds.

11) Once the restoration is seated, each cement surface/margin may be light cured for 40 seconds or allowed to self cure for 10 minutes. Note: For porcelain and pre-cured composite restorations, each cement surface/margin must be light cured for 40 seconds.

12) Instruct patient to avoid applying any pressure for 10-15 minutes.

### Instructions for bonding endodontic posts:

1) Prepare the endodontically treated tooth to receive the post (a root apex sealer and gutta percha filling approximately one third of the root canal are recommended). Trial fit and adjust post as needed. Bond to cast posts can be enhanced by using an air abrasion system and applying RelyX Ceramic Primer. Dry for 5 seconds.

2) Etching: Apply Scotchbond etchant to the prepared tooth. Wait 15 seconds. Rinse for 10 seconds. Dry for 2 seconds. Remove excess moisture with an absorbent paper point.

3) Adhesive: Apply a uniform coat to etched enamel and dentin. Remove excess pooled adhesive with absorbent paper point. Air thin for five seconds to evaporate solvents.

4) Light cure for 10 seconds. (A light transmitting post of appropriate size may be used for additional curing).

5) Dispense the appropriate amount of cement onto a mixing pad and mix for 10 seconds.

6) Apply cement to the bonding surface of the preparation (in and around canal using a periodontal probe). Place a thin layer of mixed cement on post.

7) Seat the post. While holding in place remove excess cement. Light cure for 40 seconds from the occlusal surface to allow immediate placement of core build-up material

### Instructions for bonding amalgam to both structures:

1) Isolation: Rubber dam is the recommended method of isolation.

2) Cavity preparation: Prepare a standard amalgam cavity preparation. Roughen residual restorative materials with an air abrasion system or a bur.

3) Matrix application: Lightly lubricate the inner surface of the matrix band with hard wax or petroleum jelly before placement.

4) Etching: Apply Scotchbond etchant to enamel, dentin and any residual restorative. Wait for 15 seconds. Rinse for 10 seconds. Blot excess water using a cotton pellet or mini-sponge. The surface should appear glistening without pooling of water.

5) Adhesive: Immediately after blotting, apply 2-3 consecutive coats of adhesive to etched enamel, dentin and any residual restorative material for 15 seconds with gentle agitation using a fully saturated applicator. Gently air thin for five seconds to evaporate solvents.

6) Light cure for 10 seconds.

7) Dispense the appropriate amount of cement onto a mixing pad and mix for 10 seconds.

8) Use a brush or appropriate applicator to place cement in adhesive-sealed preparation. Triturate amalgam during placement of cement.

9) Condense and burnish amalgam in the usual way.

10) Instruct patient to avoid applying any pressure for 10-15 minutes.

### Additional notes:

1) Brush handles can be disinfected in the same manner as nonimmersible handpieces, air/water syringe and ultrasonic scalers following American Dental Association (ADA) and Center for Disease Control (CDC) recommendations. (Council on Dental Materials, Instruments and Equipment and Council on Dental Therapeutics. Infection control recommendations for the dental office and the dental laboratory. JADA 116(2):241-248, 1988.)

### Storage and use:

1) Adper Scotchbond 1 XT system can be stored at room temperature.

2) Adper Scotchbond 1 XT adhesive should be capped immediately after use to minimize evaporation.

3) Do not expose materials to elevated temperature or intense light.

4) Do not store products in proximity to eugenol containing products.

5) This system is designed to be used at room temperature of approximately 21-24°C or 70-75°F.

6) Shelf life of the unit dose is 24 months at room temperature. Shelf life of the vial delivery is 36 months at room temperature. See outer package for expiry date.

7) For cleanup, Scotchbond etchant can be removed with water, while the uncured adhesive can be removed with alcohol.

No person is authorized to provide any information which deviates from the information provided in this instruction sheet.

### Warranty

3M ESPE warrants this product will be free from defects in material and manufacture. 3M ESPE MAKES NO OTHER WARRANTIES INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and 3M ESPE's sole obligation shall be repair or replacement of the 3M ESPE product.

### Limitation of Liability

Except where prohibited by law, 3M ESPE will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

### 44-0007-4225-2-B

Dental Products  
3M ESPE  
9/05

Care is required with any bonding agent used beneath precision castings because added film thickness may preclude accurate seating. Avoid adhesive pooling in areas of the preparation that would effect the fit of any prosthetic device.

Instructions for bonding crowns, bridges (including resin-bonded bridges), inlays and onlays:

- Remove temporary restoration. Trial-fit the final restoration with light finger pressure to evaluate the fit, shade and marginal integrity. Adjust if necessary.
- Prepare the bonding surface of the indirect restoration and the core build up, if applicable.

Porcelain bonding surfaces should have been etched with hydrofluoric acid by the dental laboratory. Metal and amalgam bonding surfaces should be roughened, preferably using an air abrasion system, diamond or bur. Any composite surfaces should be roughened with a diamond, bur or air abrasion system. Glass ionomer build-ups should be primed with a slurry of plain flour of pumice.

- Silane treatment (porcelain or porcelain/metal indirect restorations): Apply RelyX Ceramic Primer to the bonding surface of the indirect restoration. Dry for 5 seconds.
- Clear the prepared teeth in preparation for seating and bonding using a plain flour of pumice slurry. Rinse and dry thoroughly, isolate from moisture and adjacent teeth.

5) Etching: Apply Scotchbond etchant to both enamel and dentin. Wait 15 seconds. Rinse for 10 seconds. Blot excess water using a cotton pellet or mini-sponge. The surface should appear glistening without pooling of water.

6) Adhesive: Immediately after blotting, apply 2-3 consecutive coats of adhesive to etched enamel and dentin for 15 seconds with gentle agitation using a fully saturated applicator. Gently air thin for five seconds to evaporate solvents being careful to avoid excess adhesive on all prepared surfaces.

7) Light cure preparation for 10 seconds per bonding surface.

8) Dispense the appropriate amount of cement onto a mixing pad and mix for 10 seconds.

9) Apply and evenly distribute a thin layer of cement to the bonding surface of the indirect restoration.

10) Slowly seat and hold restoration in proper occlusion. Begin clean-up of excess cement approximately 3-5 minutes after seating. Optional: If excess cement is removed immediately after seating, each cement surface/margin must be light cured for 40 seconds.

11) Once the restoration is seated, each cement surface/margin may be light cured for 40 seconds or allowed to self cure for 10 minutes. Note: For porcelain and pre-cured composite restorations, each cement surface/margin must be light cured for 40 seconds.

12) Instruct patient to avoid applying any pressure for 10-15 minutes.

Instructions for bonding endodontic posts:

- Prepare the endodontically treated tooth to receive the post (a root apex sealer and gutta percha filling approximately one third of the root canal are recommended). Trial fit and adjust post as needed. Bond to cast posts can be enhanced by using an air abrasion system and applying RelyX Ceramic Primer. Dry for 5 seconds.

2) Etching: Apply Scotchbond etchant to the prepared tooth. Wait 15 seconds. Rinse for 10 seconds. Dry for 2 seconds. Remove excess moisture with an absorbent paper point.

3) Adhesive: Apply a uniform coat to etched enamel and dentin. Remove excess pooled adhesive with absorbent paper point. Air thin for five seconds to evaporate solvents.

4) Light cure for 10 seconds. (A light transmitting post of appropriate size may be used for additional curing).

5) Dispense the appropriate amount of cement onto a mixing pad and mix for 10 seconds.

6) Apply cement to the bonding surface of the preparation (in and around canal using a periodontal probe). Place a thin layer of mixed cement on post.

7) Seat the post. While holding in place remove excess cement. Light cure for 40 seconds from the occlusal surface to allow immediate placement of core build-up material

8) Dispense the appropriate amount of cement onto a mixing pad and mix for 10 seconds.

9) Apply and evenly distribute a thin layer of cement to the bonding surface of the indirect restoration.

10) Slowly seat and hold restoration in proper occlusion. Begin clean-up of excess cement approximately 3-5 minutes after seating. Optional: If excess cement is removed immediately after seating, each cement surface/margin must be light cured for 40 seconds.

11) Once the restoration is seated, each cement surface/margin may be light cured for 40 seconds or allowed to self cure for 10 minutes. Note: For porcelain and pre-cured composite restorations, each cement surface/margin must be light cured for 40 seconds.

12) Instruct patient to avoid applying any pressure for 10-15 minutes.

Instructions for bonding endodontic posts:

- Prepare the endodontically treated tooth to receive the post (a root apex sealer and gutta percha filling approximately one third of the root canal are recommended). Trial fit and adjust post as needed. Bond to cast posts can be enhanced by using an air abrasion system and applying RelyX Ceramic Primer. Dry for 5 seconds.

2) Etching: Apply Scotchbond etchant to the prepared tooth. Wait 15 seconds. Rinse for 10 seconds. Dry for 2 seconds. Remove excess moisture with an absorbent paper point.

3) Adhesive: Apply a uniform coat to etched enamel and dentin. Remove excess pooled adhesive with absorbent paper point. Air thin for five seconds to evaporate solvents.

4) Light cure for 10 seconds. (A light transmitting post of appropriate size may be used for additional curing).

5) Dispense the appropriate amount of cement onto a mixing pad and mix for 10 seconds.

6) Apply cement to the bonding surface of the preparation (in and around canal using a periodontal probe). Place a thin layer of mixed cement on post.

7) Seat the post. While holding in place remove excess cement. Light cure for 40 seconds from the occlusal surface to allow immediate placement of core build-up material

8) Dispense the appropriate amount of cement onto a mixing pad and mix for 10 seconds.

9) Apply and evenly distribute a thin layer of cement to the bonding surface of the indirect restoration.

10) Slowly seat and hold restoration in proper occlusion. Begin clean-up of excess cement approximately 3-5 minutes after seating. Optional: If excess cement is removed immediately after seating, each cement surface/margin must be light cured for 40 seconds.

11) Once the restoration is seated, each cement surface/margin may be light cured for 40 seconds or allowed to self cure for 10 minutes. Note: For porcelain and pre-cured composite restorations, each cement surface/margin must be light cured for 40 seconds.

12) Instruct patient to avoid applying any pressure for 10-15 minutes.

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ENGLISH –

**General Information:**  
Adper™ Scotchbond™ 1 XT Adhesive, manufactured by 3M ESPE, is a simple, moist bonding adhesive containing 10%, 5mm colloidal filler.

Adper Scotchbond 1 XT adhesive offers the dental practitioner a wide range of applications. These include bonding to all classes of direct composite restorations as well as procedures involving porcelain, composite, metal repair, et al amalgam, root surface desensitization and bonding of porcelain veneers with RelyX™ Veneer Cement System and RelyX™ Ceramic Primer, manufactured by 3M ESPE.

After light curing under Adper Scotchbond 1 XT adhesive, it may also be used for amalgam and indirect bonding procedures when combined with RelyX™ ARC adhesive resin cement, manufactured by 3M ESPE. Compatibility with indirect bonding procedures is due to the low film thickness (approximately 10µm) of cured Adper Scotchbond 1 XT.

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**Use of etchant is critical for both enamel and dentin surfaces.**

**Recommendations:**  
Use Vitrebond™ Light Cure Glass Ionomer Liner/Base, manufactured by 3M ESPE, in areas of deep cavity excavation such as Class I and II restorations. If pulp exposure has occurred, use a minimum amount of calcium hydroxide followed by an application of Vitrebond liner/base. Adper Scotchbond 1 XT adhesive will bond to Vitrebond liner/base whether or not the ionomer was treated with etchant.

Adper Scotchbond 1 XT adhesive includes an etch of enamel and dentin as a part of the procedure. It is recommended that the surfaces be left moist during rinsing. Excess surface moisture should be removed by blotting.

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#### Precautions for Dental Personnel and Patients

**Adper Scotchbond™ Etchant, manufactured by 3M ESPE, contains 35% weight w/v phosphoric acid.**

Protective eyewear for patients and dental staff is recommended when using etchants. Avoid contact with oral soft tissue, eyes and skin. If accidental eye contact occurs, flush immediately with large amounts of water. For eye contact, consult a physician.

**Adper Scotchbond 1 XT adhesive contains acrylates including HEMA (2-hydroxyethylmethacrylate).** Avoid use of this product on patients with known acrylic allergies. To reduce the risk of allergic response, minimize exposure to these materials. In particular, avoid exposure to uncured resins. Use of **protective gloves and a new wax technique is recommended.** If skin contact occurs, wash skin with soap and water. Acrylates may penetrate commonly-used gloves. If adhesive contacts glove, remove and discard glove, wash hands immediately with soap and water and then re-glove. If accidental contact with eyes or prolonged contact with oral soft tissue occurs, flush with large amounts of water. If irritation persists, consult a physician. See Material Safety Data Sheet (MSDS) for additional cautionary information. You may obtain the current MSDS by visiting the website: <http://www.3m.com/MSDS> or contacting your 3M ESPE Dental Products representative.

**Sensitivity:**  
Some patients may experience transitory postoperative sensitivity. The risk of sensitivity can be minimized by the following measures:

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Adhesive Application:

Use of compressed air is not recommended to remove pooled water remaining after the etch step–Blot excess moisture from the preparation using a cotton pellet or mini-sponge.

Adhesive immediately after blotting.

**Restorative:**  
Place restorative material in increments, curing each increment separately.

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Adjust occlusion carefully. Check for hyperocclusion, particularly in lateral excusion contacts.

#### Etchant Syringe Assembly:

- 1) Protective eyewear for patients and staff is recommended when using the delivery system.
- 2) Prepare the delivery system: Remove cap from etchant syringe and SAVE. Twist a blue disposable tip securely onto the syringe. Holding the syringe with the tip away from the patient and any dental staff, express a small amount of etchant onto a dispensing pad or a 2 x 2 gauze to assure that the delivery system is not clogged. If clogged, remove the dispensing tip and express a small amount of etchant directly from the syringe. Remove any visible plug, if present, from the syringe opening. Replace dispensing tip and again express etchant. If plug remains, discard dispensing tip and replace with a new one. Bore the dispensing tip to a desired angle. Place bend midway along tip. Do not bend dispensing tip at its hub as this may cause the cap to break free.
- 3) Delivery system storage: Remove used dispensing tip and discard. Twist on storage cap. **Storage of the delivery syringe with a used dispensing tip or without storage cap will allow drying of the etchant and consequent clogging of the system.** Replace storage cap with a new dispensing tip at next use.
- 4) If desired, the etchant may be extruded onto a dispensing pad and applied with a brush or other appropriate instrument.
- 5) If a liquid etchant is desired, the etchant may be dispensed into a dappen dish and stirred to increase its fluidity.
- 6) Disinfection: Discard used dispensing tip. Replace syringe cap. Disinfect the capped syringe in the same manner as nonimmisable handpieces, air/water syringe and ultrasonic scalers following American Dental Association (ADA) and Centers for Disease Control (CDC) recommendations. (Council on Dental Materials, Instruments and Equipment and Council on Dental Therapeutics. Infection control recommendations for the dental office and the dental laboratory. JADA 116(2):241-248, 1988.)

#### Dispensing Adper Scotchbond 1 XT Adhesive

**Unit Dose:**  
Attention: To minimize risk of accidental contact with eyes and skin, hold the shaft of the disposable applicator over the opening where the applicator enters the foil package with your thumb and index finger. Do not activate the unit dose without a disposable applicator. With your thumb and index finger of the other hand, squeeze the large blister to transfer the adhesive into the channel enclosing the applicator. Briefly spin the applicator to fully saturate with adhesive.

Vial:

Pinch the sides of the cap to release the locking mechanism and flip the cap back to reveal the dispensing tip. Squeeze the exact number of drops you need into the disposable mixing well. When finished, flip the cap back until it is secured by its locking mechanism.

#### Instructions for direct light cure restorations in enamel and dentin:

- 1) Isolation: Rubber dam is the preferred method of isolation.
- 2) Cavity preparation: Prepare cavity with minimal tooth reduction. Bevel cavosurface enamel margins.
- 3) Etching: Apply Scotchbond etchant to both enamel and dentin. Wait 15 seconds. Rinse for 10 seconds. Blot excess water using a cotton pellet or mini-sponge. The surface should appear glistening without pooling of water.
- 4) Adhesive: Immediately after blotting, apply 2-3 consecutive coats of adhesive to etched enamel and dentin for 15 seconds with gentle agitation using a fully saturated applicator. Gently air thin for five seconds to evaporate solvents. Light cure for 10 seconds.
- 5) Restorative placement, cure and finishing: Refer to manufacturer's instructions for placement, cure and finishing of restorative materials.

#### Instructions for bonding porcelain veneers:

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- Do for 5 seconds.
- 2) Clean the prepared teeth in preparation for seating and bonding using a plain four of pumice slurry. Rinse thoroughly and dry.
- 3) Try in veneer with RelyX™ Try-In Paste, manufactured by 3M ESPE. After try in, isolate from adjacent teeth with clear matrix strip.
- 4) Etching: Apply Scotchbond etchant to both enamel and dentin. Wait 15 seconds. Rinse for 10 seconds. Blot excess water using a cotton pellet or mini-sponge. The surface should appear glistening without pooling of water. An additional 15 second etch may be appropriate for teeth that were not prepared using a diamond or bur. Residual organic matter can also require additional etch time.
- 5) Adhesive: Immediately after blotting, apply 2-3 consecutive coats of adhesive to etched enamel and dentin for 15 seconds with gentle agitation using a fully saturated applicator. Gently air thin for five seconds to evaporate solvents. Do not light cure.
- 6) Adhesive application to veneer: Apply 1 coat of adhesive to the acid etched, silane treated veneer. Dry thoroughly. Do not light cure.
- 7) Luting material application to veneer: Apply RelyX Veneer Cement to the bonding surface of veneer.
- 8) Seating and curing: Carefully seat the veneer. Clean excess luting cement from the veneer margins. Cure each area of the veneer for times recommended by the luting cement manufacturer. We recommend curing the gingival margin first, followed by the body and the incisal margin. Avoid direct contact with the light-guide during initial curing.

#### Instructions for bonding to composite and self amalgam:

- 1) Isolation: Rubber dam is the preferred method of isolation.
- 2) Roughen the existing material: Roughen the surface of existing amalgam a composite using either a bur, diamond or a sandblast technique.

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- 4) Adhesive: Immediately after blotting, apply 2-3 consecutive coats of adhesive to etched enamel, dentin and existing restorative material for 15 seconds with gentle agitation using a fully saturated applicator. Gently air thin for five seconds to evaporate solvents. Light cure for 10 seconds.
- 5) Masking: In the case of self amalgam, mask the metal surface with a thin liner (0.25 - 0.5mm) of the appropriate 3M ESPE Masking Agent paste, manufactured by 3M ESPE, with a brush. Light cure each layer for 20 seconds.
- 6) Restorative placement: Refer to manufacturer's instructions for placement, cure and finishing of restorative material.

#### Instructions for porcelain repair:

- 1) Isolation: Rubber dam is the preferred method of isolation.
- 2) Preparation: Clean the surface to be repaired with a slurry of plain four of pumice. Rinse and dry thoroughly. Roughen the surface of existing metal or porcelain using either a bur, diamond or a sandblast technique. Be careful to remove all loose porcelain and bevel the margin. Remove surface glaze using a plain four of pumice and the margin.
- 3) Etching: Apply Scotchbond etchant to all substrates. Wait for 15 seconds. Rinse for 10 seconds. Dry 5 seconds.
- 4) Silane treatment for porcelain and metal: Apply RelyX Ceramic Primer (No. 2721) to the etched surface and dry.
- 5) Adhesive: Apply 2 consecutive coats of Adper Scotchbond 1 XT adhesive to silane treated porcelain or metal. Dry gently for 5 seconds. Light cure for 10 seconds.
- 6) Masking: To opacity the metal before the final composite placement, mask the metal surface with a thin layer (0.25 - 0.5mm) of the appropriate 3M ESPE Masking Agent paste using a brush. Light cure each layer for 20 seconds.
- 7) Restorative placement, cure and finishing: Refer to manufacturer's instructions for placement, cure and finishing of restorative material.

#### Instructions for root surface desensitization:

- 1) Lightly rough the root surface with four of pumice. Rinse and blot dry.
- 2) Etching: Apply Scotchbond etchant to the etched dentin. Wait for 15 seconds. Rinse for 10 seconds. Blot excess water using a cotton pellet or mini-sponge. The surface should appear glistening without pooling of water.
- 3) Adhesive: Immediately after blotting, apply 2-3 consecutive coats of adhesive to etched enamel and dentin for 15 seconds with gentle agitation using a fully saturated applicator. Gently air thin for five seconds to evaporate solvents. Light cure for 10 seconds. Apply 2 additional coats of adhesive. Dry gently for 5 seconds. Light cure for 10 seconds.
- 4) Remove the oxygen inhibited layer with a moistened gauze.

#### Instructions for Amalgam and Indirect Bonding Procedures Using RelyX ARC Adhesive resin cement:

Physical properties of today's esthetic indirect restorations require that they be bonded into place to maximize the strength of the restoration as well as the tooth. A general perception may exist that light cure adhesives cannot be used for indirect restorations. It's true that many conventional light cure adhesives have a higher film thickness and can not be used under a food protected However, Adper Scotchbond 1 XT adhesive is ethanol/water based, has a low film thickness (approximately 10µm) and should not interfere with the seating of indirect restorations.

#### Note:

Care is required with any bonding agent used beneath precision castings because added film thickness may preclude accurate seating. **Adap, adhesive pooling in areas of the preparation that would effect the fit of any prosthetic device.**

#### Instructions for bonding crowns, bridges (including resin-bonded bridges), implants and onlays:

- 1) Remove temporary restoration. Trim/fit the final restoration with light finger pressure to evaluate the fit, shade and marginal integrity. Adjust if necessary.
- 2) Prepare the bonding surface of the indirect restoration and the core build up, if applicable. Porcelain bonding surfaces should have been etched with hydrofluoric acid by the dental laboratory. Metal and amalgam bonding surfaces should be roughened, preferably using an air abrasion system, diamond or bur. Any composite surfaces should be roughened with a diamond, bur or air abrasion system. Glass ionomer bulk cements should be pumiced with a slurry of plain four of pumice.
- 3) Silane treatment (porcelain or porcelain/metal indirect restorations): Apply RelyX Ceramic Primer to the bonding surface of the indirect restoration. Dry for 5 seconds.
- 4) Clean the prepared teeth in preparation for seating and bonding using a plain four of pumice slurry. Rinse and dry thoroughly. Isolate from moisture and adjacent teeth.
- 5) Etching: Apply Scotchbond etchant to both enamel and dentin. Wait 15 seconds. Rinse for 10 seconds. Blot excess water using a cotton pellet or mini-sponge. The surface should appear glistening without pooling of water.
- 6) Adhesive: Immediately after blotting, apply 2-3 consecutive coats of adhesive to etched enamel and dentin for 15 seconds with gentle agitation using a fully saturated applicator. Gently air thin for five seconds to evaporate solvents being careful to avoid excess adhesive on all prepared surfaces.
- 7) Light cure preparation for 10 seconds per bonding surface.
- 8) Dispense the appropriate amount of cement onto a mixing pad and mix for 10 seconds.
- 9) Apply and evenly distribute a thin layer of cement to the bonding surface of the indirect restoration.
- 10) Slowly seat and hold restoration in proper occlusion. Begin clean-up of excess cement approximately 3-5 minutes after seating. **Optional: If excess cement is removed immediately after seating, each cement surface/margin must be light cured for 40 seconds.**

- 11) Once the restoration is seated, each cement surface/margin may be light cured for 40 seconds or allowed to self cure for 10 minutes. **Note: for porcelain and pre-cured composite restorations, each cement surface/margin must be light cured for 40 seconds.**
  - 12) Instruct patient to avoid applying any pressure for 10-15 minutes.
- Instructions for bonding endosteic posts:**
- 1) Prepare the endosteically treated tooth to receive the post (a root apex sealer and gutta percha filling approximately one third of the root canal are recommended). **Trim the fit and adjust post as needed.** Bond to cast posts can be enhanced by using an air abrasion system and then applying RelyX Ceramic Primer. Dry for 5 seconds.
  - 2) Etching: Apply Scotchbond etchant to the prepared tooth. Wait 15 seconds. Rinse for 10 seconds. Dry for 2 seconds. Remove excess moisture with an absorbent paper point.
  - 3) Adhesive: Apply a uniform coat to etched enamel and dentin. Remove excess pooled adhesive with absorbent paper point. Air thin for five seconds to evaporate solvents.
  - 4) Light cure for 10 seconds. (A light transmitting spot of appropriate size may be used for additional curing.)
  - 5) Dispense the appropriate amount of cement onto a mixing pad and mix for 10 seconds.
  - 6) Apply cement to the bonding surface of the preparation (in and around canal using a periodontal probe). Place a thin layer of mixed cement on post.
  - 7) Seat the post. While holding in place remove excess cement. Light cure for 40 seconds from the occlusal surface to allow immediate placement of core build-up material

#### Instructions for bonding amalgam to tooth structure:

- 1) Isolation: Rubber dam is the recommended method of isolation.
  - 2) Cavity preparation: Prepare a standard amalgam cavity preparation. Roughen residual restorative materials with an air abrasion system or a bur.
  - 3) Matrix application: Lightly lubricate the inner surface of the matrix band with hard wax or petroleum jelly before placement.
  - 4) Etching: Apply Scotchbond etchant to enamel, dentin and any residual restorative. Wait for 15 seconds. Rinse for 10 seconds. Blot excess water using a cotton pellet or mini-sponge. The surface should appear glistening without pooling of water.
  - 5) Adhesive: Immediately after blotting, apply 2-3 consecutive coats of adhesive to etched enamel, dentin and any residual restorative material for 15 seconds with gentle agitation using a fully saturated applicator. Gently air thin for five seconds to evaporate solvents.
  - 6) Light cure for 10 seconds.
  - 7) Dispense the appropriate amount of cement onto a mixing pad and mix for 10 seconds.
  - 8) Use a brush or appropriate applicator to place cement in adhesive-sealed preparation.
  - Trabete amalgam during placement of cement.**
  - 9) Condense and burnish amalgam in the usual way.
  - 10) Instruct patient to avoid applying any pressure for 10-15 minutes.
- Additional notes:**
- 1) Brush handles can be disinfected in the same manner as nonimmisable handpieces, air/water syringe and ultrasonic scalers following American Dental Association (ADA) and Center for Disease Control (CDC) recommendations. (Council on Dental Materials, Instruments and Equipment and Council on Dental Therapeutics. Infection control recommendations for the dental office and the dental laboratory. JADA 116(2):241-248, 1988.)

#### Storage and use:

- 1) Adper Scotchbond 1 XT system can be stored at room temperature.
- 2) Adper Scotchbond 1 XT adhesive should be capped immediately after use to minimize evaporation.
- 3) Do not expose materials to elevated temperature or intense light.
- 4) Do not store products in proximity to eugenol containing products.
- 5) This system is designed to be used at room temperature of approximately 21-24°C or 70-75°F.
- 6) Shelf life of the unit dose is 24 months at room temperature. Shelf life of the vial delivery is 36 months at room temperature. See outer package for expiry date.
- 7) For cleanup, Scotchbond etchant can be removed with water, while the uncured adhesive can be removed with alcohol.

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#### Warranty

3M ESPE warrants this product will be free from defects in material and manufacture. 3M ESPE MAKES NO OTHER WARRANTIES INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. User is responsible for determining the suitability of the product for use's application. If this product is defective within the warranty period, your exclusive remedy and 3M ESPE's sole obligation shall be repair or replacement of the 3M ESPE product.

**Limitation of Liability**  
Except where prohibited by law, 3M ESPE will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

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