


**Flowable** when desired.  
**Stable** as required.



## The Tetric EvoCeram Family of Composites

Tetric EvoFlow is a component of the Tetric line of products. Whether you are in need of a sculptable or a flowable composite, this line offers you a suitable option for any indication.

SCULPTABLE

Traditional  
2mm layering

Bulk filling

**Tetric EvoCeram®**  
*The composite*

- 22 Universal shades
- The composite that has proven longevity and over 90 million restorations world wide.

**Tetric EvoCeram® Bulk Fill**  
*The bulk-fill composite*



- 3 Posterior Shades
- 4-mm bulk-fill composite for posterior restorations

FLOWABLE

Incremental layering

**Tetric EvoFlow®**  
*The flowable composite*

- 14 Shades 3 new shades for all flowable indications
- The time-tested flowable composite with high radiopacity



Pictures: Dr Eduardo Mahn, Santiago, Chile

## New shades to complete the system


The shade range has been extended to include three new shades: B1 Enamel, A2 Dentin and B2 Dentin. Shade B1 is widely used for lighter teeth while the dentin materials can be used to mask light stains in posterior teeth.

Tetric EvoFlow is ideal for use as a liner in combination with other composites such as Tetric EvoCeram, Tetric EvoCeram Bulk Fill or any other regular-viscosity resin based composites.



## Exceptional radiopacity

High radiopacity is a prerequisite for a sound diagnosis of secondary caries. Tetric EvoFlow exhibits a radiopacity of 360% AI and makes it easier for clinicians to properly distinguish between dental tissue and restorative material.



R&D Ivoclar Vivadent AG, Schaan, 2012 (measured acc. to ISO 4049) \*These trademarks are not registered trademarks of Ivoclar Vivadent AG.

## Both flowable and stable

When lining a restoration, restoring a Class V or restoring another small restoration, a flowable that remains stable and stays in place is ideal. Tetric EvoFlow is characterized by optimum surface affinity and simultaneous stability. It flows where it is needed and stable when required.



Tetric EvoFlow (left), SDR Flow\* (middle) and Palfique Estelite LV\* (right) have been applied to a vertically positioned glass slab.

R&D Ivoclar Vivadent AG, Schaan, 2012.

\*This trademark is not a registered trademark of Ivoclar Vivadent AG.

## Ideal consistency to suit all indications



Liner under regular-viscosity composite restorations Class I to V, (e.g. Tetric EvoCeram/Tetric EvoCeram Bulk Fill restorations)



Preventive resin restorations (sealant)



Small restorations (Class III and V)



## Delivery Forms



## KITS

**595951 Tetric EvoFlow Syringe Assortment Kit**  
 Contents:  
 8 - syringes (2 g each) Tetric EvoFlow  
 (1 each in shades; A1, A2, A3, A3.5, A4, B3,  
 A3.5 Dentin, Translucent)  
 Shade Guide

**638381 Tetric EvoFlow Cavifil Assortment Kit**  
 Contents:  
 40 - cavifils (0.2 g each) Tetric EvoFlow  
 (5 each in shades; A1, A2, A3, A3.5, A4, B3,  
 A3.5 Dentin, Translucent)  
 Shade Guide

## REFILLS

### Syringe Refills (1 x 2g)

595953 Tetric EvoFlow Refill 1x2 g A1  
 595954 Tetric EvoFlow Refill 1x2 g A2  
 595955 Tetric EvoFlow Refill 1x2 g A3  
 595956 Tetric EvoFlow Refill 1x2 g A3.5  
 595978 Tetric EvoFlow Refill 1x2 g A4  
 641358 Tetric EvoFlow Refill 1x2 g B1  
 595957 Tetric EvoFlow Refill 1x2 g B3  
 595979 Tetric EvoFlow Refill 1x2 g T  
 641360 Tetric EvoFlow Refill 1x2 g A2 Dentin  
 595980 Tetric EvoFlow Refill 1x2 g A3.5 Dentin  
 641361 Tetric EvoFlow Refill 1x2 g B2 Dentin  
 595981 Tetric EvoFlow Refill 1x2 g Bleach XL  
 595982 Tetric EvoFlow Refill 1x2 g Bleach L  
 595984 Tetric EvoFlow Refill 1x2 g Bleach M

### Cavifil Refills (20 x 0.2g)

595987 Tetric EvoFlow Refill 20x0.2 g A1  
 595988 Tetric EvoFlow Refill 20x0.2 g A2  
 595989 Tetric EvoFlow Refill 20x0.2 g A3  
 595990 Tetric EvoFlow Refill 20x0.2 g A3.5  
 595992 Tetric EvoFlow Refill 20x0.2 g A4  
 641348 Tetric EvoFlow Refill 20x0.2 g B1  
 595991 Tetric EvoFlow Refill 20x0.2 g B3  
 595993 Tetric EvoFlow Refill 20x0.2 g T  
 641349 Tetric EvoFlow Refill 20x0.2 g A2 Dentin  
 595994 Tetric EvoFlow Refill 20x0.2 g A3.5 Dentin  
 641350 Tetric EvoFlow Refill 20x0.2 g B2 Dentin  
 595995 Tetric EvoFlow Refill 20x0.2 g Bleach XL  
 595996 Tetric EvoFlow Refill 20x0.2 g Bleach L  
 596008 Tetric EvoFlow Refill 20x0.2 g Bleach M

## SHADE GUIDE

**595952** Tetric EvoFlow Shade Guide

## LUERLOCK TIPS

**634849** Luerlock Tips 0.9-mm Black Qty. 20

100% CUSTOMER SATISFACTION  
 GUARANTEED!  
[ivoclarvivadent.com](http://ivoclarvivadent.com)

©2012 Ivoclar Vivadent, Inc.  
 Tetric EvoFlow and Tetric EvoCeram are  
 registered trademarks of Ivoclar Vivadent.  
 646476 Rev 1 9/2012

**United States**  
 175 Pineview Drive  
 Amherst, New York 14228  
 800-533-6825  
 Fax: 716-691-2285  
[www.ivoclarvivadent.com](http://www.ivoclarvivadent.com)

**Canada**  
 1 - 6600 Dixie Road  
 Mississauga, Ontario L5T 2Y2  
 800-263-8182  
 Fax: 905-670-3102  
[www.ivoclarvivadent.com](http://www.ivoclarvivadent.com)

**ivoclar  
vivadent®**  
 passion vision innovation