

Your practice is *our* inspiration.™

TempBond | Strong, reliable, short-term attachment



TempBond

Clinically proven performance you trust



TempBond™ temporary cement for trial restorations or temporary crowns and bridges withstands mastication, seals the restoration and prevents seepage, yet is easy to remove. Non-eugenol TempBond NE™ will not inhibit the polymerization of permanent resin cements and acrylic temporaries. And resin-based TempBond Clear™ with Triclosan is the most translucent provisional cement in the market for superior esthetics.

Inspired by:

"In this age of cosmetic dentistry, TempBond Clear is the temporary cement for provisionally luting anterior temporary restorations into place."

– Dr. Gary Radz

Research

TempBond was developed in 1961 based on zinc oxide-eugenol chemistry, a classic setting reaction with over 100 years of clinical history. Then in 1989, Kerr introduced TempBond NE with eugenol-free zinc oxide, substituting carboxylic acids in place of eugenol. And TempBond Clear with Triclosan incorporates a unique redox initiator system and low refractive index fillers for excellent color stability and high translucency.



Case courtesy of Dr. Tony Soileau

Features

Delivery choices. Unidose®, tubes and new automix syringe offer flexibility in delivery options.

Flows and mixes easily. Delivers optimal consistency for solid, complete seating of restorations.

High bond strength. Heightens patient confidence. Prevents leakage and sensitivity to cold and heat.

Easy to remove. Enhances patient comfort. Minimizes potential for damage to the prep or temporary.



Ordering

- 00370 TempBond Original Tube**
Contains: 1 Tube Base (50 g), 1 Tube Accelerator (15 g) and 1 Mixing Pad
- 31377 TempBond Original Unidose**
Contains: 50 Unidose Packets (2.4 g each) and 1 Mixing Pad
- 33215 TempBond Original Automix Syringe**
Contains: 2 Syringes (11.8 g each) and 20 Mixing Tips
- 21370 TempBond NE Tube**
Contains: 1 Tube Base (50 g), 1 Tube Accelerator (15 g) and 1 Mixing Pad
- 31973 TempBond NE Unidose**
Contains: 50 Unidose Packets (2.4 g each) and 1 Mixing Pad
- 33217 TempBond NE Automix Syringe**
Contains: 2 Syringes (11.7 g each) and 20 Mixing Tips
- 33351 TempBond Clear with Triclosan Automix Syringe**
Contains: 1 Syringe (6 g) and 10 Mixing Tips

Accessories

- 28366 Mixing Pad (60 Sheets – 2 3/8" x 1 7/8")
- 33361 Automix Tip Refill (50 Tips)

Associated Products

Fill-In™ Temporary Crown and Bridge Material

- 32496 Fill-In Intro Kit
Contains: 5 Syringes (3 g each) of Shades: A1, A2, A3.5, XL and Instruction Card

Fill-In Refill Kits

Each Refill Kit contains 35 syringes (3 g each)

- 32491 A1 32493 – A3.5
- 32492 A2 32494 – XL

Tempphase™ Temporary Crown and Bridge Material Refills

Each Refill contains 1 Cartridge (72 g) and 12 Mix Tips

Regular Set Fast Set

- 29366 A1 29367 – A1
- 28070 A2 28655 – A2
- 28071 A3.5 28656 – A3.5
- 28072 B1 28657 – B1
- 28073 C2 28658 – C2

Tempphase Accessories

- 28747 Extruder (1) 28075 Mix Tips (48)

©2008 Kerr Corporation Lit. No. 34123