




Ultradent Products, Inc.

# TISSUE MANAGEMENT

— Since 1978 —



**Rapid, profound hemostasis**  
*is paramount for proper tissue management.*


With a complete line of tissue management solutions, Ultradent enables you to control bleeding and sulcular fluid, an imperative step in creating quality direct and indirect restorations.



Bleeding prep




Syringe with brush tip



Illustrated brush tip with capillaries



Non-bleeding/clean prep



## VISCOSTAT® CLEAR

25% Aluminum Chloride

- For the esthetic zone, no staining
- Viscous—won't run, drip, or migrate to other teeth

*Great for same-day placement of restoration*



## VISCOSTAT®

20% Ferric Sulfate

- Associated with routine use
- For rapid, profound hemostasis
- Viscous—won't run, drip, or migrate to other teeth

*Excellent as an alternative to formocresol for vital pulpotomies*



## ASTRINGEDENT®

15.5% Ferric Sulfate

- Stops moderate bleeding in seconds
- Also eliminates sulcular fluid contamination for optimal bonding

*The "Classic" hemostatic agent—Ultradent's first product*



## ASTRINGEDENT® X

12.7% Iron Solution

Facilitates profound hemostasis, even with the most challenging cases, when used with Ultradent's Metal Dento-Infusor tip.

- Controls very problematic, difficult-to-stop bleeding
- Less acidic than competitive iron subsulfate solutions

*For especially difficult cases of bleeding*



Metal Dento-Infusor Tip



Blue Mini Dento-Infusor Tip

## DENTO-INFUSOR TIPS

For clean, dry preparations

- Essential to rapid, profound hemostasis and sulcular fluid control
- Padded brushes infuse hemostatic agents into bleeding capillaries and wipe coagulum away

*Dento-Infusor tips offer direct precision and placement*