


simple, adaptable



endodontic solutions

# NaviTip Delivery Tips

When it comes to root canals, Ultradent is with you every step of the way



NaviTip® FX® Tip

NaviTip® Tip

NaviTip® Double Sideport

**Don't change your technique. Make it easier with NaviTip.**

Ultradent's NaviTip tips are the optimal delivery tips for every step in a root canal procedure. Unlike traditional irrigation tips, the cannula of each NaviTip is slightly rigid through the base and center, but flexible at the tip to allow for easy navigation of curved canals.

## Simplify any procedure

Use NaviTip tips to deliver any endodontic irrigant directly where it's needed—just short of the apex.




**ULTRADENT**  
PRODUCTS, INC.

800.552.5512 [ULTRADENT.COM](http://ULTRADENT.COM)



### NaviTip® 29ga or 30ga

- Controlled delivery anywhere in the canal
- Flexible cannula easily navigates curved canals



### NaviTip® 31ga Double Sideport Irrigator Tip

- The small cannula navigates the most intricate canals
- Double sideports horizontally express irrigant toward the canal walls, not the apex



### NaviTip® 30ga FX® Tips

- Brush tip scrubs, cleans, and irrigates simultaneously
- Flexible cannula easily navigates curved canals
- Listed as an "EXCELLENT" product by a prominent independent research institute.\*

### NaviTip® Tips

	20pk		50pk	
	29ga	30ga	29ga	30ga
17mm, White	#5112	#1249	#1378	#1421
21mm, Yellow	#5113	#1349	#1374	#1422
25mm, Blue	#5114	#1250	#1376	#1423
27mm, Green	#5115	#1354	#1377	#1424
Assorted	#5116	#1351	#1379	#3319

### NaviTip® Double Sideport Irrigator Tip

	20pk	50pk
21mm, Yellow	#5121	#5122
27mm, Green	#5123	#5124

### NaviTip® FX® Tips

	20pk
17mm, White	#1452
25mm, Blue	#1454

## NaviTip® adapts to your style

The NaviTip tip delivers low-viscosity materials exactly where they are needed—even in the smallest canals.

Thanks to NaviTip 31ga's double sideport openings, your irrigant is always safely expressed against the canal walls rather than toward, or even through, the foramen.



NaviTip® Tip



NaviTip®  
Double Sideport

## Use NaviTip to deliver products listed below and many other products:

- **ChlorCid®** & **ChlorCid® V** Sodium Hypochlorite Solutions
- **Ultradent® EDTA** 18% Solution
- **Consepsis®** & **Consepsis® V** Chlorhexidine Antibacterial Solutions
- **File-Eze®** EDTA Lubricant
- **Ultradent® Citric Acid** 20% Solution
- **UltraCal® XS** Calcium Hydroxide Paste
- **EndoREZ®** Hydrophilic Resin Root Canal Sealer



\*Clinical Research Associates Newsletter. Volume 29, Issue 1, January 2005.